**IFU**

Product: EURODERM INCISE IODINE

INDICATION FOR USE

Euroderm Incise Iodine is intended to be used as mechanical barrier for surgical incisions. It contains Iodine which grants a continuous broad-spectrum antimicrobial activity, through the surgical procedure.

CONTRA-INDICATIONS

Do not use in case of known intolerance to the material of support (PU) or to the Iodine.

DESCRIPTION OF PRODUCT

Euroderm Incise Iodine it is a sterile, adhesive, transparent, antistatic, non-glare incision drape in polyurethane film which makes clearly visible the incision site.

The film, permeable to gas and moisture, but impermeable to liquids and bacteria, protects the area from any contamination and allows breathe, preventing moisture accumulate under the drape.

Euroderm Incise Iodine securely adheres up to the wound edge for all the time requested by the operation and, furthermore, it properly holds the drapes around the incision area, reducing so the use of other fixing devices.

The elasticity, flexibility and conformability of the film make it perfect for easy applications, even around difficult anatomical contours.

The film is resistant but can be cleanly cut with the scalpel.

Euroderm Incise Iodine is latex free.

The device,

according to the MDD 93/42 EEC of 14.06.1993 (Annex IX, Rule 13), belongs to class III Sterile.

PACKAGING

Each box, labelled according to MDD 93/42 EEC, contains 10 pcs. individually pouched.

The stability of the packaging has been tested as conform to the criteria of UNI EN 868 and grants to the medical devices a validity of sterility of 5 years since the date of sterilization.

REF. SIZES AND PACKAGING:

Article	Sizes	Adhesive area	Packaging
717910	cm. 10 x cm. 24	cm. 10 x cm. 18	10 pcs ind.pouched
717915	cm. 15 x cm. 26	cm. 15 x cm. 20	10 pcs ind.pouched
717930	cm. 30 x cm. 26	cm. 30 x cm. 20	10 pcs ind.pouched
717940	cm. 40 x cm. 34	cm. 40 x cm. 28	10 pcs ind.pouched
717942	cm. 42 x cm. 42	cm. 42 x cm. 36	10 pcs ind.pouched
717944	cm. 44 x cm. 26	cm. 44 x cm. 20	10 pcs ind.pouched
717945	cm. 45 x cm. 55	cm. 45 x cm. 49	10 pcs ind.pouched
717950	cm. 50 x cm. 90	cm. 50 x cm. 84	10 pcs ind.pouched
717960	cm. 60 x cm. 52	cm. 60 x cm. 46	10 pcs ind.pouched
717961	cm. 60 x cm. 35	cm. 60 x cm. 29	10 pcs ind.pouched
717970	cm. 70 x cm. 80	cm. 70 x cm. 74	10 pcs ind.pouched
717980	cm. 80 x cm. 60	cm. 80 x cm. 54	10 pcs ind.pouched
717990	cm. 90 x cm. 80	cm. 90 x cm. 74	10 pcs ind.pouched
717991	cm. 90 x cm. 50	cm. 90 x cm. 44	10 pcs ind.pouched
717992	cm. 90 x cm. 90	cm. 90 x cm. 84	10 pcs ind.pouched

IDENTIFICATION OF RISK:

Injurious effects on health: not determinable

Injurious effects on the environment: Does not represent risks for the environment provided that the national laws for the refuses are respected.

Ways of fire extinguishing:

Suitable: Foam, Powder, and CO₂.

Contra-indicated: none of our knowledge, next to fire suitable ways of fire extinguishing have to be used.

Special packaging conditions: None.

Rev. 01 – 05/2022

WARNING:

Single use
Do not use if the individual package is damaged or opened.
Do not use the product after the expiration date.

INSTRUCTIONS FOR USE:

1. Wash and disinfect hands thoroughly.
2. Disinfect the surgical site according to the hospital policy.
3. Open the pouch.
4. Unwind the drape and, keeping it in tension from both edges, remove the protective paper of the adhesive side (2-3 persons depending on the drape size).
5. Apply the drape on the incision site avoiding wrinkles.
6. Remove the protective paper under the blue stripes in order to make firmly adhere the sides of the drape.

MANUFACTURING OF THE PRODUCT:

Euroderm Incise Iodine is manufactured in a clean room (Class 100.000) under the following conditions:

Temperature: 18 – 25 °C

Relative humidity: 30% - 80%

BIOBURDEN: The control is currently made on one sample every 4 batches according to UNI EN 11737-1:2006 and the microbe burden must be inferior to 100 CFU/sample.

PROCESS OF STERILIZATION:

The article is sterilised with ethylene oxide (10% EO – 90% CO₂) in our works. The entire process of sterilisation, from the preconditioning to the degassing, has been confirmed and monitored in accordance to the Norm UNI EN ISO 1135-1:2008.

STORAGE:

Keep dry and far from heating sources.
Do not damage the packaging in order ensure its hygienic security.
Shelf life: 5 years

The re-validation is carried out every three years.

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