




Pain Numbing Applicator



Makes perfect sense to lose that sense of pain

Numbs site prior to any type of injection procedure

- An alternative to dermal creams
- Highly cost effective per use
- Simple to operate
- Reusable
- Over 50,000 Hospitals & Doctors are using Coolsense worldwide

Paediatrics
Emergency
Operating theatres
Day Surgery
Oncology
Dialysis
Needle-phobic patients



Coolsense

Coolsense is a unique device that numbs the site prior to injection and can be used in place of dermal anesthetic creams or ice. It is an excellent option to use prior to all types of intramuscular and Intravenous injections such as PRP, cortisone injections, blood collection/pathology, local and general anesthesia, routine immunization, hormone injections, dialysis, chemotherapy, and much more.

- **A reusable device** that is applied to the skin for up to 3 seconds (face) or a maximum of 10 seconds (body), after which time you can administer an injection/cannulate immediately
- **Simple and quick to apply** without the mess and lengthy time required for dermal anesthetic creams such as EMLA to work. Instead of applying cream to multiple sites you just apply the device to one site at a time without causing vascular constriction